

Sulfobutylether Beta-Cyclodextrin as Chiral Selector in Capillary Electrophoresis

Cyclodextrins (CDs) are the most frequently used chiral selectors in the capillary electrophoresis [1]. For academic research usually single isomer CDs are used, but for practical purposes the random substituted CDs give often suitable results especially when CD derivatives produced industrially are used. Sulfobutylether beta-CD (SBE- β -CD) is a pharmaceutical excipient in US Pharmacopoeia as solubilizing agent [2] and is a component of several drug formulations in the market. Its composition is strictly regulated therefore the average degree of substitution (DS) and the distribution of the components of various DS fall always in a narrow range. The standard quality results in high reproducibility for analytical applications.

The first report on the application of SBE- β -CD in capillary electrophoresis (CE) was published in 1994 by the group of Valentino Stella [3]. In this pioneering work, the authors demonstrated that enhanced enantiomeric separation can be achieved when the electrophoretic mobility of the chiral selector is opposite to that of the analyte. The "extended" enantioselectivity ability of SBE-CDs, based on the advantage of the countercurrent flow of the negatively charged additive with respect to the electroosmotic flow (EOF), was then further deepened by other scientists. Along the years, the family of sulfobutylated CDs gained popularity and SBE- β -CD is often referred as the most versatile chiral additive among the negatively charged CDs (and CDs in general) as it allows the enantioseparation of a wide variety of pharmacologically active racemates.

The use of SBE-CDs in CE is suitable for the analysis of illicit amphetamine, methamphetamine, methcathinone and propoxyphene [4] and for the enantiomeric separation of ephedrine and related compounds [5]. Chankvetadze et al. described the use of SBE- β -CD as a chiral additive for the resolution of basic racemic drugs, such as clenbuterol, dimethindene, etilefrine, mefloquine and metomidate, in free capillary zone electrophoresis.



The resolution of the racemates was achieved with very low concentrations (micromolar) of the chiral additive and the high efficiency of the selector was attributed to its counter-current mobility in respect to the racemic solute [6].

Also Desiderio and Fanali reported on the effective use of SBE- β -CD as a chiral selector in CE in a study on the enantiomeric separation of warfarin, pindolol, propranolol, terbutaline, etc., a variety of underivatized anionic and cationic compounds of pharmaceutical interest as well as dansyl-amino acids [7].

SBE- β -CD was found as best chiral selector for the enantiomeric separation of several neutral and anionic herbicides by capillary zone electrophoresis [7]. SBE- γ -CD has been also investigated as chiral selector for CE. A comparative study between SBE- γ -CD and native γ -CD has shown that SBE- γ -CD is more efficient in the electrophoretic separation of cationic analytes, but it cannot generally replace γ -CD since the enantioselectivities of the two selectors are sometimes distinctly different [8]. Electrophoretic chiral separation of chloroquine and pemoline was effectively achieved by using 2.5 mM SBE- β -CD and 1.0 mM SBE- β -CD in 50 mM pH 2.5 sodium phosphate buffer, respectively [9]. Fundamental drugs used for the clinical treatment of Parkinson's disease, such as DOPA and structurally related compounds were resolved into their enantiomers by CE with SBE- β -CD as chiral selector [10–11]. The use of SBE- β -CD as chiral additive in CE was also useful for the separation of estrogens and amino acid derivatives [12–14]. The enantioseparations of the antimalarian erythro-mefloquine and its analogues were screened by CE using library of CD derivatives, including SBE- β -CD, as chiral selectors [15].

The chiral separation of the dihydropyridine calcium channel blockers (DHPs) with CE has been reported utilizing SBE- β -CD as the chiral selectors: it was found that the anionic SBE- β -CD gives a robust separation of amlodipine enantiomers when using acetonitrile - 20 mM NaH_2PO_4 (pH 3.95), containing 20 mM SBE- β -CD (35:100, v/v) [16]. Kong et al. developed an electrophoretic method for the enantioseparation of neutral m-nisoldipine, by comparing SBE- β -CD and carboxymethyl- β -CD (CM- β -CD) as chiral selectors. It was concluded that SBE- β -CD was superior chiral selector than CM- β -CD. It was also pointed out that the DS of the SBE- β -CD can remarkably influence the enantiorecognition process. In particular, high DS SBE- β -CD induced better investigated by CD-mediated CE and it was found that the enantiomers of the neutral dihydropyridines were baseline-separated only with SBE- β -CD [18].

Fillet et al. performed comparative chiral resolution studies among the commonly used neutral CDs with SBE- β -CD. Increased enantiomeric resolution for various acidic drugs was achieved by using dual systems comprising neutral β -CD and charged SBE- β -CD [19].

The extensive work on the use of anionic CDs in chiral CE have been summarized in excellent reviews [20–22]. De Boer et al. in their review reported several examples about the utility of sulfobutylated β - and γ -CDs as chiral selectors and the authors referred to SBE- β -CD



as the "most widely used anionic CD derivative" [23]. In the review by Mikuš et al. [24], some new achievements in the use of SBE-CDs are reported.

Further papers on the analytical application of SBE-CDs are summarized in Table 1.

Table 1. Examples of using SBE- β -CD as chiral selector

Analytes	References
tadalafil ^a	25
sitagliptin ^b	26
imperanene ^c	27
R-modafinil	28
ofloxacin and its metabolites	29
volinanserin and its intermediate	30
methylargininies	31
polycyclic aromatic hydrocarbons	32-34
bupivacaine hydrochloride enantiomers	35
isradipine enantiomers	36
chloroquine enantiomers	37
benzofurys, cathinones, diphenidines, ethylphenidate, methiopropamine and thiothinone	38
brompheniramine, chlorpheniramine, cetirizine and promethazine racemate ^d	39
pantoprazole enantiomers	40
retroamide enantiomers	41
alogliptin enantiomers	42, 43
trans-tramadol and its metabolites ^e	44, 45
pyroglutamic acid derivatives ^e	46

^a SBE- α -CD, ^b SBE- β -CD (average DS \sim 4), ^c SBE- γ -CD (DS \sim 4), ^d SBE- β -CD and native β -CD dual CD system, ^e electrokinetic chromatography (EKC)

The versatility of SBE-CDs as capillary electrophoretic tools, together with their commercial availability largely justify the extensive research conducted so far on this family of derivatives. The large number and wide structural variety of guest molecules where efficient chiral resolution was achieved demonstrate not only the excellent enantio-recognition but also the general and high affinity complex forming ability of these cyclodextrin derivatives.



References

1. L. Szente, J. Szemán, Cyclodextrins in analytical chemistry: host-guest type molecular recognition, *Anal. Chem.*, 85 (2013) 8024–8030
2. I. Puskás, E. Varga, K. Tuza, J. Szemán, É. Fenyvesi, T. Sohajda, L. Szente, Sulfobutylether-cyclodextrins: structure, degree of substitution and functional performance, In: *Cyclodextrins*, Chapter 10, ISBN: 978-1-63482-788-1, Editor: F.G. Ramirez © 2015 Nova Science Publishers, Inc.
3. R. Tait, D. O. Thompson, V. Stella, J. F. Stobaugh, Sulfobutyl ether β -cyclodextrin as a chiral discriminator for use with capillary electrophoresis. *Anal. Chem.* 66 (1994) 4013-4018
4. I. S. Lurie, R. F. X. Klein, T. A. Dal Cason et al., Chiral resolution of cationic drugs of forensic interest capillary electrophoresis with mixtures of neutral and anionic cyclodextrins. *Anal. Chem.* 66 (1994) 4019–4026
5. C. Dette, S. Ebel, S. Terabe, Neutral and anionic cyclodextrins in capillary zone electrophoresis: enantiomeric separation of ephedrine and related compounds, *Electrophoresis* 15 (1994) 799–803
6. B. Chankvetadze, G. Endresz, G. Blaschke, About some aspects of the use of charged cyclodextrins for capillary electrophoresis enantioseparation. *Electrophoresis* 15(6) (1994) 807
7. C. Desiderio, S. Fanali, Use of negatively charged sulfobutyl ether- β -cyclodextrin for enantiomeric separation by capillary electrophoresis, *J. Chromatogr. A* 716 (1995) 183-196
8. M. Jung, E. Francotte, Comparison of γ -cyclodextrin sulfobutyl ether and unmodified γ -cyclodextrin as chiral selectors in capillary electrophoresis, *J. Chromatogr. A* 755 (1996) 81–88
9. Y. Dong, A. Huang, Y. Sun, Z. Sun, Chiral separation of chloroquine and pemoline by capillary zone electrophoresis with sulfobutyl ether β -cyclodextrin as buffer additive, *Chromatographia* 48 (1998) 310-313
10. M. Dolezalová, S. Fanali, Enantiomeric separation of dihydroxyphenylalanine (DOPA), methyl-dihydroxyphenylalanine (MDOPA) and hydrazinomethyl-dihydroxyphenylalanine (CDOPA) by using capillary electrophoresis with sulfobutyl ether- β -cyclodextrin as a chiral selector, *Electrophoresis* 21 (2000) 3264-3269
11. S. Sarac, B. Chankvetadze, G. Blaschke, Enantioseparation of 3,4-dihydroxyphenylalanine and 2-hydrazino-2-methyl-3-(3,4-dihydroxyphenyl)propanoic acid by capillary electrophoresis using cyclodextrins. *J. Chromatogr. A* 875 (2000) 379–387
12. Y. Deng, M. Huang, Capillary electrophoretic separation and theoretical study of inclusion complexes of sulfobutyl ether β -cyclodextrin with estrogens, *Int. J. Quantum Chem.* 100 (2004) 746-752
13. Y. Deng, J. Zhou, S. M. Lunte et al, Capillary electrophoretic separation of estrogens using anionic sulfobutyl ether β -cyclodextrin, *Anal. Commun.* 34 (1997) 129-131
14. C. Perrin, M.G. Vargas, Y. Vander Heyden, M. Maftouh, D. L. Massart, Fast development of separation methods for the chiral analysis of amino acid derivatives using capillary electrophoresis and experimental designs, *J. Chromatogr. A* 883 (2000) 249 -265
15. B. Chankvetadze, N. Burjanadze, G. Blaschke, Enantioseparation of erythromefloquine and its analogues in capillary electrophoresis, *J. Pharm. Biomed. Anal.* 32 (2003) 41-49
16. P. K. Owens, A. F. Fell, M. W. Coleman et al., Effect of charged and uncharged chiral additives on the resolution of amlodipine enantiomers in liquid chromatography and capillary electrophoresis. *J. Chromatogr. A* 797 (1998) 187–195
17. D. Kong, W. Yang, K. Duan, Y. Cui, S. Xing, X. Zhang, X. Sheng, L. Zhang, Capillary electrophoretic enantioseparation of m-nisoldipine using two different β cyclodextrins, *J. Sep. Sci.* 32 (2009) 3178–3183
18. T. Cristians, U. Holzgrabe, Enantioseparation of dihydropyridine derivatives by means of neutral and negatively charged β -cyclodextrin derivatives using capillary electrophoresis, *Electrophoresis* 21 (2000) 3609–3717
19. M. Fillet, I. Bechet, J. Crommen et al., Enantiometric separation of acidic drugs by capillary electrophoresis using a combination of charged and uncharged β -cyclodextrins as chiral selectors, *J. High Resol. Chromatogr.* 19 (1996) 669–673
20. B. Chankvetadze, G. Endresz, G. Blaschke, Charged cyclodextrin derivatives as chiral selectors in capillary electrophoresis, *Chem. Soc. Rev.* 25 (1996) 141–153
21. J. M. Saz, M. L. Marina, Recent advances in the use of cyclodextrin in the chiral analysis of drugs by capillary electrophoresis, *J. Chromatogr. A* 1467 (2016) 79–94
22. W. Tang, S.-C. Ng, D. Sun (Eds.), *Modified cyclodextrin for chiral separation*, Springer, 2013 paragraph: 6.4 Sulfobutyl Ether- β -CD (SBE- β -CD), pp 209-216
23. T. De Boer, R. A. De Zeeuw, G. J. De Jong, K. Esing, Recent innovations in the use of charged cyclodextrins in capillary electrophoresis for chiral separations in pharmaceutical analysis,



- Electrophoresis 21 (2000) 3220–3239
24. P. Mikuš, K. Marákova, Advanced CE for chiral analysis of drugs, metabolites, and biomarkers in biological samples, *Electrophoresis* 30 (2009) 2773–2802
 25. I. Fejős, A. Kazsoki, T. Sohajda, E. Márványos, B. Volk, L. Szente, S. Béni, Interactions of non-charged tadalafil stereoisomers with cyclodextrins: Capillary electrophoresis and nuclear magnetic resonance studies, *J. Chromatogr. A*, 1363 (2014) 348–55
 26. W. Wang, S. Xiang, B. Xiang et al., Enantiomeric separation and determination of the enantiomeric impurity of armodafinil by capillary electrophoresis with sulfobutyl ether- β -cyclodextrin as chiral selector, *Molecules* 17 (2012) 303–314
 27. T. Sohajda, Z. Szakacs, L. Szente, B. Noszal, S. Beni, Chiral recognition of imperanene enantiomers by various cyclodextrins: a capillary electrophoresis and NMR spectroscopy study, *Electrophoresis* 33 (2012) 1458–1464
 28. T. Sohajda, W. Hu, S. Béni et al., Evaluation of the interaction between sitagliptin and cyclodextrin derivatives by capillary electrophoresis and nuclear magnetic resonance spectroscopy, *Electrophoresis* 32 (2011), 2648–2654
 29. C. Horstkotter, G. Blaschke, Stereoselective determination of ofloxacin and its metabolites in human urine by capillary electrophoresis using laser-induced fluorescence detection, *J. Chromatogr. B* 754, (2001) 169–178
 30. L. Escuder-Gilabert, Y. Martìn-Biosca, M. J. Medina-Hernàndez, S. Sagrado, Cyclodextrins in capillary electrophoresis: Recent developments and new trends, *J. Chromatogr. A* 1357 (2014) 2–23
 31. T. H. Linz, S. M. Lunte, Heat-assisted extraction for the determination of methylarginines in serum by CE, *Electrophoresis* 34 (2013) 1693–1700
 32. L. Ferey, N. Delaunay, D. N. Rutledge, A. Huertas, Y. Raoul, P. Gareil, J. Vial, Use of response surface methodology to optimize the simultaneous separation of eight polycyclic aromatic hydrocarbons by capillary zone electrophoresis with laser-induced fluorescence detection, *J. Chromatogr. A* 1302 (2013) 181–190
 33. L. Ferey, N. Delaunay, D. N. Rutledge, C. B. Y. Cordella, H. This, A. Huertas, Y. Raoul, P. Gareil, Optimizing separation conditions of 19 polycyclic aromatic hydrocarbons by cyclodextrin-modified capillary electrophoresis and applications to edible oils, *Talanta* 119 (2014) 572–581
 34. L. Ferey, N. Delaunay, D. N. Rutledge, A. Huertas, Y. Raoul, P. Gareil, J. Vial, I. Rivals, An experimental design based strategy to optimize a capillary electrophoresis method for the separation of 19 polycyclic aromatic hydrocarbons, *Anal. Chim. Acta* 820 (2014) 195–204
 35. Y. Li, Y. J. Yu, P. Y. Zhu, G. L. Duan, Y. Li, F. Y. Song, Chiral separation of bupivacaine hydrochloride by capillary electrophoresis with high frequency conductivity detection and its application to rabbit serum and pharmaceutical injection, *Pharmazie* 67 (2012) 25–30
 36. F. A. Aguiar, C. M. de Gaitani, K. B. Borges, Capillary electrophoresis method for the determination of isradipine enantiomers: Stability studies and pharmaceutical formulation analysis, *Electrophoresis* 32 (2011) 2673–2682
 37. S. Wongwan, G. K. E. Scriba, Development and validation of a capillary electrophoresis assay for the determination of the stereoisomeric purity of chloroquine enantiomers, *Electrophoresis* 32 (2011) 2669–2672
 38. Taschwer M, Hofer MG, Schmid MG, Enantioseparation of benzofurys and other novel psychoactive compounds by CE and sulfobutylether β -cyclodextrin as chiral selector added to the BGE, *Electrophoresis* 35 (2014) 2793–2799
 39. Z. I. Szabó, C. Tóth, G. Hancu, D. L. Muntean, Simultaneous chiral separation of four H1 antihistamines by capillary zone electrophoresis using a dual cyclodextrin system, *Chromatographia* 78 (2015) 1377–1384
 40. J. Guan, F. Yan, S. Shi, S. Wang, Optimization and validation of a new CE method for the determination of pantoprazole enantiomers, *Electrophoresis* 33 (2012) 1631–1636
 41. D. Baudelet, C. Furman, A. Ghinet, X. Dezitter, S. Adriouch, F. Capet, T. Rogez-Florent, P. Gautret, B. Rigo, R. Millet, C. Vaccher, E. Lipka, Evaluation and comparison of three different separation techniques for analysis of retroamide enantiomers and their biological evaluation against h-P2X7 receptor, *J. Chromatogr. B* 986–987 (2015) 35–43
 42. I. Fejős, Z. Urbancsok, W. Zhou, T. Sohajda, W. Hu, L. Szente, S. Béni, Separation of alogliptin enantiomers in cyclodextrin-modified capillary electrophoresis: a validated method, *Electrophoresis* 35 (2014) 2885–2891
 43. A. Kazsoki, I. Fejős, T. Sohajda, W. Zhou, W. Hu, L. Szente, S. Béni, Development and validation of a cyclodextrin-modified capillary electrophoresis method for the enantiomeric separation of vildagliptin enantiomers, *Electrophoresis* 37 (2016) 1318–1325
 44. S. Rudaz, S. Cherkaoui, P. Dayer, S. Fanali, J. L. Veuthey, Simultaneous stereoselective analysis of tramadol and its main phase I metabolites by on-line capillary zone electrophoresis-electrospray



ionization mass spectrometry, *J. Chromatogr. A* 868 (2000) 295–303

45. H. C. Liu, Y. Yu, N. Wang, M. Deng, J. F. Liu, H. Y. Xue, Pharmacokinetics of the enantiomers of trans-tramadol and its active metabolite, trans-o-demethyltramadol, in healthy male and female Chinese volunteers, *Chirality* 16 (2004) 112–118
46. D. Baudelet, A. Ghinet, C. Furmann, X. Dezitter, P. Gautret, B. Rigo, R. Millet, C. Vaccer, E. Lipka, Antagonists of the P2X7 receptor: mechanism of enantioselective recognition using highly sulfated and sulfobutylether cyclodextrins by capillary electrokinetic chromatography, *Electrophoresis* 35 (2014) 2892–2899

Milo Malanga & Éva Fenyvesi

CycloLab Cyclodextrin R&D Laboratory, Ltd.,
Budapest, HUNGARY



Bibliography & Keywords of Selected Publications of the Month

Arima, H.; Motoyama, K.; Higashi, T.

Potential Use of Cyclodextrins as Drug Carriers and Active Pharmaceutical Ingredients

Review

Chemical & pharmaceutical bulletin, 2017, 65, 341-348;

Arima, H.; Motoyama, K.; Higashi, T.

Potential therapeutic application of dendrimer/cyclodextrin conjugates with targeting ligands as advanced carriers for gene and oligonucleotide drugs

Polyamidoamine dendrimer (G3) conjugates with α -CyD, Cell-specific drug carriers, Polyethylene glycol, galactose, lactose, mannose, fucose and folic acid-appended α -CDEs

Therapeutic Delivery, 2017, 8, 215-232; DOI:10.4155/tde-2016-0064

Shi, X.; Li, W.; Liu, H.; Yin, D.; Zhao, J.

β -Cyclodextrin induces the differentiation of resident cardiac stem cells to cardiomyocytes through autophagy

Regeneration of damaged hearts, Increased expression of cardiac transcription factors and structural proteins, Enhanced myogenic transition, Cholesterol efflux

Biochimica et Biophysica Acta, Molecular Cell Research, 2017, 1864, 1425-1434; DOI:10.1016/j.bbamcr.2017.05.012

Perez-Anes, A.; Szarpak-Jankowska, A.; Jary, D.; Auzely-Velty, R.

β -CD-Functionalized Microdevice for Rapid Capture and Release of Bacteria

Quartz-crystal microbalance, Functionalization of surfaces by β -CD, Selectively methylated β -CD derivative as a competitive host, Capture/elution strategy, Escherichia coli

ACS Applied Materials & Interfaces, 2017, 9, 13928-13938; DOI:10.1021/acsami.7b02194

Zhai, Z.; Zhang, F.; Chen, X.; Zhong, J.; Liu, G.; Tian, Y.; Huang, Q.

Uptake of silver nanoparticles by DHA-treated cancer cells examined by surface-enhanced Raman spectroscopy in a microfluidic chip

Cyclodextrin-coated silver nanoparticles modified with para-aminothiophenol and folic acid, Endocytosis, Targeted drug delivery, Controlled release

Lab on a Chip, 2017, 17, 1306-1313; DOI:10.1039/C7LC00053G

Sun, H-L.; Chen, Y.; Han, X.; Liu, Y.

Tunable Supramolecular Assembly and Photoswitchable Conversion of Cyclodextrin/Diphenylalanine-Based 1D and 2D Nanostructures

Supramolecular assembly of adamantanyl-modified diphenylalanine with azobenzene-bridged bis(β -cyclodextrin), Nanosheet, Nanotube, Photo-switchable one-dimensional/two-dimensional morphological interconversion

Angewandte Chemie, International Edition, 2017, Ahead of Print; DOI:10.1002/anie.201612629



Fine-Shamir, N.; Beig, A.; Zur, M.; Lindley, D.; Miller, J.; M. Dahan, A.

Toward Successful Cyclodextrin Based Solubility-Enabling Formulations for Oral Delivery of Lipophilic Drugs: Solubility-Permeability Trade-Off, Biorelevant Dissolution, and the Unstirred Water Layer

Success/failure of cyclodextrin-based solubility-enabling formulations for oral delivery of low-solubility drugs, Danazol, HP β CD, Excess CD

Molecular Pharmaceutics, 2017, Ahead of Print; DOI:10.1021/acs.molpharmaceut.7b00275

Oster, M. and Schlatter, G. and Gallet, S. and Baati, R. and Pollet, E. and Gaillard, C. and Averous, L. and Fajolles, C. and Hebraud, A.

The study of the pseudo-polyrotaxane architecture as a route for mild surface functionalization by click chemistry of poly(ϵ -caprolactone)-based electrospun fibers

Pseudo-polyrotaxanes of cyclodextrin and PCL, Core:shell fibers, Coaxial electrospinning, Grafting fluorescein isothiocyanate and bicyclononyne groups

Journal of Materials Chemistry B: Materials for Biology and Medicine, 2017, 5, 2181-2189; DOI:10.1039/C6TB03089K

Liu, S.; Zhong, C.; Chen, J.; Zhan, J.; He, J.; Zhu, Y.; Wang, Y.; Wang, L.; Ren, L.

Thermoresponsive Self-Assembled β -Cyclodextrin-Modified Surface for Blood Purification

Removal of bilirubin from the body in clinics, Detoxification system, Adamantane on gold surface, β -CD dimer, Resistance to plasma proteins, In vitro platelet adhesion assay, Hemolysis assay, Negligible cytotoxicity

ACS Biomaterials Science & Engineering, 2017, Ahead of Print; DOI:10.1021/acsbmaterials.7b00156

Benkovics, G.; Malanga, M.; Fenyvesi, E.

The 'Visualized' macrocycles: Chemistry and application of fluorophore tagged cyclodextrins

Review, chemosensing, Fluorescent microscopy, Synthetic strategies, Photodynamic therapy

International Journal of Pharmaceutics, 2017, Ahead of Print; DOI:10.1016/j.ijpharm.2017.04.035

Terauchi, M.; Inada, T.; Kanemaru, T.; Ikeda, G.; Tonegawa, A.; Nishida, K.; Arisaka, Y.; Tamura, A.; Yamaguchi, S.; Yui, N.

Potentiating bioactivity of BMP-2 by polyelectrolyte complexation with sulfonated polyrotaxanes to induce rapid bone regeneration in a mouse calvarial defect

Sulfopropyl ether (SPE)-modified α -cyclodextrins threaded along a poly(ethylene glycol) chain

Journal of Biomedical Materials Research, Part A, 2017, 105, 1355-1363; DOI:10.1002/jbm.a.36016

Singh, V.; Guo, T.; Wu, L.; Xu, J.; Liu, B.; Gref, R.; Zhang, J.

Template-directed synthesis of a cubic cyclodextrin polymer with aligned channels and enhanced drug payload

CD-cubes, metal organic frameworks, γ -CDs in MOFs crosslinked by diphenyl carbonate, Enhanced adsorption capacity for doxorubicin

RSC Advances, 2017, 7, 20789-20794; DOI:10.1039/C7RA02306E

Kitagishi, H.; Minegishi, S.

Iron(II)porphyrin-Cyclodextrin Supramolecular Complex as a Carbon Monoxide-Depleting Agent in Living Organisms

Per-O-methylated β -cyclodextrin (TMe- β -CD), HemoCD

Chemical & pharmaceutical bulletin, 2017, 65, 336-340



Motoyama, K.; Nishiyama, R.; Maeda, Y.; Higashi, T.; Ishitsuka, Y.; Kondo, Y.; Irie, T.; Era, T.; Arima, H.

Synthesis of multi-lactose-appended β -cyclodextrin and its cholesterol-lowering effects in Niemann-Pick type C disease-like HepG2 cells

Active targeting hepatocytes with β -CD, Decreasing the concentration of intracellular cholesterol, Negligible cytotoxicity, Increased internalization

Beilstein Journal of Organic Chemistry, 2017, 13, 10-18; DOI:10.3762/bjoc.13.2

Canterini, S.; Dragotto, J.; Dardis, A.; Zampieri, S.; Egle De Stefano, M.; Mangia, F.; Erickson, R. P.; Fiorenza, M. T.

Shortened primary cilium length and dysregulated Sonic hedgehog signaling in Niemann-Pick C1 Disease

Dysregulation of Shh signaling, 2-Hydroxypropyl- β -cyclodextrin, Abnormal morphogenesis of the cerebellum of Npc1-deficient mice, Primary cilium

Human molecular genetics, 2017,

Collins, C. J.; Loren, B. P.; Mondjinou, Y.; Skulsky, J. L.; Chaplain, C. R.; Alam, M. S.; Haldar, K.; Thompson, D. H.

Pluronic based β -cyclodextrin polyrotaxanes for treatment of Niemann-Pick Type C disease

HP- β -CD pro-drug, polyrotaxanes

Scientific reports, 2017, 7, 46737

Massaro, M.; Colletti, C. G.; Lazzara, G.; Guernelli, S.; Noto, R.; Riela, S.

Synthesis and Characterization of Halloysite-Cyclodextrin Nanosponges for Enhanced Dyes Adsorption

Inorganic-organic nanosponge hybrids, Clay, Microwave irradiations in solvent-free conditions, Rhodamine B adsorption

ACS Sustainable Chemistry & Engineering, 2017, 5, 3346-3352; DOI:10.1021/acssuschemeng.6b03191

Rey-Rico, A.; Babicz, H.; Cucchiari, M.; Madry, H.; Concheiro, A.; Alvarez-Lorenzo, C.

Supramolecular polypseudorotaxane gels for controlled delivery of rAAV vectors in human mesenchymal stem cells for regenerative medicine

Cartilage regeneration, Hyaluronic acid, Chondroitin sulfate, α CD, Improved transgene expression. Sustained release

International journal of pharmaceuticals, 2017

Nakahata, M.; Takashima, Y.; Harada, A.

Supramolecular Polymeric Materials Containing Cyclodextrins

Review, Macroscopic self-assembly of polymer gels carrying host and guest residues, Stimuli-responsive self-healing properties, Macroscopic motion of artificial muscle

Chemical & pharmaceutical bulletin, 2017, 65, 330-335

Letort, S.; Bosco, M.; Cornelio, B.; Bregier, F.; Gouhier, G.; Estour, F.; Daulon, S.

Structure-efficiency relationships of cyclodextrin scavengers in the hydrolytic degradation of organophosphorus compounds

Organophosphorus nerve agents, Competitive inclusion between the organophosphorus substrate and the iodosobenzoate group, Imidazole group

Beilstein journal of organic chemistry, 2017, 13, 417-427



Gogolashvili, A.; Tatumashvili, E.; Chankvetadze, L.; Sohajda, T.; Szeman, J.; Salgado, A.; Chankvetadze, B.

Separation of enilconazole enantiomers in capillary electrophoresis with cyclodextrin-type chiral selectors and investigation of structure of selector-selectand complexes by using nuclear magnetic resonance spectroscopy

Enantiomer migration order, Heptakis(2-O-methyl-3,6-di-O-sulfo)- β -CD, Shallow external complex

Electrophoresis, 2017, Ahead of Print; DOI:10.1002/elps.201700078

Patyk-Kazmierczak, E.; Warren, M. R.; Allan, D. R.; Katrusiak, A.

Pressure inverse solubility and polymorphism of an edible γ -cyclodextrin-based metal-organic framework

Trigonal structure

Physical Chemistry Chemical Physics, 2017, 19, 9086-9091; DOI:10.1039/C7CP00593H

Martel, B.; Campagne, C.; Behary, M. N.

When textiles help your recovery

Wound dressings, Implantable devices, Advanced spinning technologies of biostable or bioresorbable polymers and surface treatment technologies, Cyclodextrins releasing the active agents in situ

Medicine sciences : M/S, 2017, 33, 73-80

Castriciano, M. A.; Zagami, R.; Casaletto, M. P.; Martel, B.; Trapani, M.; Romeo, A.; Villari, V.; Sciortino, M. T.; Grasso, L.; Guglielmino, S.; Scolaro, L. M.; Mazzaglia, A.

Poly(carboxylic acid)-Cyclodextrin/Anionic Porphyrin Finished Fabrics as Photosensitizer Releasers for Antimicrobial Photodynamic Therapy

Polypropylene (PP) fabric finished with citrate-hydroxypropyl- β CD polymer, Photoinactivation of microorganisms, Staphylococcus aureus, Pseudomonas aeruginosa

Biomacromolecules, 2017, 18, 1134-1144; DOI:10.1021/acs.biomac.6b01752

Cheng, M.; Wang, M.; Zhang, S.; Liu, F.; Yang, Y.; Sun, L.; Wan, B.

Photocatalytic H₂ production using a hybrid assembly of an [FeFe]-hydrogenase model and CdSe quantum dot linked through a thiolato-functionalized cyclodextrin

Hydrogen evolution reaction, Hybrid assembly of a sulfonate-functionalized [FeFe]-hydrogenase mimic and CdSe quantum dot, 6-Mercapto- β -cyclodextrin asstabilizing reagent and a host compound for the diiron catalyst

Faraday discussions, 2017

Yuan, L-M.; Ma, W.; Xu, M.; Zhao, H-L.; Li, Y-Y.; Wang, R-L.; Duan, A-H.; Ai, P.; Chen, X-X.

Optical resolution and mechanism using enantioselective cellulose, sodium alginate and hydroxypropyl- β -cyclodextrin membranes

p-Hydroxy phenylglycine, Chiral separation, Adsorption - association - diffusion

Chirality, 2017, 29, 315-324; DOI:10.1002/chir.22693

Benkovics, G.; Afonso, D.; Sortino, S.; Darcsi, A.; Beni, S.; Conoci, S.; Fenyvesi, E.; Szente, L.; Malanga, M.

Novel β -cyclodextrin-eosin conjugates

6-Monoamino- β -cyclodextrin scaffold, Photogeneration of singlet oxygen

Beilstein journal of organic chemistry, 2017, 13, 543-551



Jicsinszky, L.; Caporaso, M.; Calcio, G. E.; Cravotto, G.; Giovannoli, C.

Synthesis of Randomly Substituted Anionic Cyclodextrins in Ball Milling

Minimizing the reagent's hydrolysis, Higher degree of substitution (DS) and more homogeneous DS distribution, Mechanochemical activation

Molecules, 2017,



Edited and produced by: CYCLOLAB

Homepage: www.cyclolab.hu

H-1525 P.O. 435, Budapest,
Hungary

Tel.: (+361)347-6060

Fax.: (+361)347-6068

e-mail: cyclolab@cyclolab.hu

