

CYCLOLAB



The Cyclodextrin Company




Characterization and application of single-isomer cyclodextrin derivatives

Gábor Benkovics


CycloLab, Cyclodextrin R&D Ltd., Budapest, Hungary

Faculty of Science, Charles University in Prague, Czech Republic


Cyclodextrins



≡




≡



CD	Number of glucose units
α	6
β	7
γ	8


Isomerism in cyclodextrin derivatization

Random/statistical substitution




-  position 2-O
-  position 3-O
-  position 6-O


Monosubstitution




Persubstitution






Disubstitution



Inclusion-assisted liquid chromatography CD-Screen columns





Inclusion-assisted ion-exchange chromatography (CD-Screen-IEC columns)



HPLC conditions: CD-Screen-IEC stationary phase, ammonium acetate (pH 4.0), gradient with ACN, ELS detection

Single-isomer carboxymethyl cyclodextrins



heptakis-(2,3-di-O-methyl-6-O-carboxymethyl)-β-cyclodextrin (HDMCM)
octakis-(2,3-di-O-methyl-6-O-carboxymethyl)-γ-cyclodextrin (ODMCM)

Preparation:




- 4-step synthesis from native CDs in 52 % (HDMCM) and 48 % (ODMCM) total yield

Characterization:


- HPLC - inclusion-assisted ion-exchange chromatography (CD-Screen-IEC columns) with ELS detection
- Capillary electrophoresis with indirect UV detection
- ¹H NMR titration

Single-isomer carboxymethyl cyclodextrins (application)

- pH tunable ionic interactions on a well-defined position




electroneutral



octaanionic


CE electropherograms applying ODMCM as chiral selector

- applicability in aqueous as well as in non-aqueous capillary electrophoresis (NACE)




NACE electropherograms applying HDMCM as chiral selector



Single-isomer sulfobutyl- β -cyclodextrin



Sulfobutylated β -CD (DS=6.4)
(Dexolve™, Captisol™)




Heptakis-(6-O-sulfobutyl)- β -CD (DS=7.0)
6-(SB)₇- β -CD





HPLC conditions: CD-Screen-IEC stationary phase, solvent gradient of 1.5% triethylamine-formate (pH 4.5) with ACN, ELS detection


Single-isomer sulfobutyl β -cyclodextrin





CE electropherograms applying 6-(SB)₇- β -CD as chiral selector





Beyond enantiomer separation: Methylated cyclodextrins




randomly-methylated β -cyclodextrins (RAMEB, CRYSMEB)



heptakis-(6-*O*-methyl)- β -cyclodextrin (SIXMEB)
heptakis-(2-*O*-methyl)- β -cyclodextrin (TWOMEB)



heptakis-(2,6-di-*O*-methyl)- β -cyclodextrin (2,6-DIMEB)
heptakis-(2,3-di-*O*-methyl)- β -cyclodextrin (2,3-DIMEB)





heptakis-(2,3,6-tri-*O*-methyl)- β -cyclodextrin (TRIMEB)



Beyond enantiomer separation



- Nanocarriers for the delivery of antimicrobials to fight resistance mechanisms
- Cyclodextrin-penicillin interactions



- ODMCM, HDMCM and TRIMEB are able to protect the β -lactam ring from catalytic degradation



Summary

- New negatively charged single-isomer CDs have been synthesized:
ODMCM , HDMCM (carboxymethylated)
6-SB₇-β-CD (sulfoalkylated)
- Chiral resolution properties of these derivatives have been studied with a wide set of test racemates applying aqueous and non-aqueous capillary electrophoresis
- The synthesized library of single-isomer carboxymethylated and methylated CD derivatives helped us to understand the inclusion phenomenon between CDs and β-lactam antibiotics
- ODMCM, HDMCM and TRIMEB are the first CD derivatives showing stabilization effect of the β-lactam ring

Acknowledgement

Milo Malanga Ph.D., Julianna Szemán Ph.D., Tamás Sohajda Ph.D.,
Éva Fenyvesi Ph.D., Erzsébet Varga, Lajos Szente Ph.D., D.Sc.
CycloLab Cyclodextrin R&D, Ltd., Budapest, Hungary

Assoc. Prof. Jindřich Jindřich
Faculty of Science, Charles University in Prague, Czech Republic

András Darcsi Ph.D., Ida Fejős Ph.D., Assoc. Prof. Szabolcs Béni
Department of Pharmacognosy, Semmelweis University, Budapest, Hungary

Agnieszka Popielec, Prof. Thorsteinn Loftsson,
Faculty of Pharmacy, University of Iceland, Iceland

European Commission, FP7 PEOPLE-ITN project n. 608407-CycloN Hit

Prof. Gyula Vígh

CYCLOLAB



The Cyclodextrin Company



Thank you for your kind attention!

**www.cyclolab.hu
benkovics@cyclolab.hu**