

MAIL: BUDAPEST, P.O.BOX 435, H-1525 HUNGARY LOCATION: ILLATOS ÚT 7. BUDAPEST, H-1097 HUNGARY TEL: +36 (1) 347 6070

E-MAIL: cyclolab@cyclolab.hu WEB: www.cyclolab.hu VAT No.: HU 10678970



Packaging Specification of Sulfobutyl-ether- β -cyclodextrin sodium salt (SBECD, Dexolvetm) bulk material

Primary Packaging Material: Polyethylene bag (sterile, gamma-irradiated)

Dimensions: thickness: 0.100 mm; 40 x 90 cm

Tests and methods	Requirements
Appearance, packaging (visual)	Colorless bag, free from visible mechanical contaminations, the packaging is free from any observable damage
Identification	The received material complies with CycloLab Purchase Order
Size (visual)	Declared dimensions ± 2 %

Storage: at room temperature; Shelf life: 5 years (the expiration date of the packaging material is determined based on the stability test)

Secondary Packaging Material: heat-sealed, three-layered foil bag two pressure relief valves

Material: inner layer: PE (60 μm), middle layer: aluminum (7 μm), outer layer: PET (12 μm)

Dimensions: thickness: 0.079 mm; 75 x 45 cm

Tests and methods	Requirements
Appearance, packaging (visual)	Free from visible mechanical contaminations, the packaging is free from any observable damage
Identification	The received material complies with CycloLab Purchase Order
Size (visual)	Declared dimensions ± 2 %

Storage: at room temperature; Shelf life: 5 years (the expiration date of the packaging material is determined based on the stability test)

CYCLO PIONEERS IN B CYCLODEXTRIN-BASED INNOVATION

CYCLOLAB CYCLODEXTRIN RESEARCH & DEVELOPMENT LABORATORY LTD.

MAIL: BUDAPEST, P.O.BOX 435, H-1525 HUNGARY LOCATION: ILLATOS ÚT 7. BUDAPEST, H-1097 HUNGARY TEL: +36 (1) 347 6070

E-MAIL: cyclolab@cyclolab.hu WEB: www.cyclolab.hu VAT No.: HU 10678970

ISO 9001

Tertiary Packaging Materials:

HDPE blue drum with locking clamp (60 L)

Material: HDPE granules; Clamp material: Steel

Dimensions: 60 x 38 cm; Total weight: 1.9 kg

Tests and methods	Requirements
Appearance (visual)	Free from visible mechanical contaminations and damages
Closure (visual)	Presence of sealing material

Storage: at room temperature.

Fiber drums with locking clamp (Small, 15 L)

Material: fiber; Clamp material: Steel

Dimensions: 31 x 20.5 cm Total weight: 2.2 kg

Tests and methods	Requirements
Appearance (visual)	Free from visible mechanical contaminations and damages
Closure (visual)	Presence of sealing material

Storage: at room temperature.

Fiber drums with locking clamp (Medium, 40 L)

Material: fiber; Clamp material: Steel

Dimensions: 40 x 34.5 cm Total weight: 3.0 kg

Tests and methods	Requirements
Appearance (visual)	Free from visible mechanical contaminations and damages
Closure (visual)	Presence of sealing material

Storage: at room temperature.



MAIL: BUDAPEST, P.O.BOX 435, H-1525 HUNGARY LOCATION: ILLATOS ÚT 7. BUDAPEST, H-1097 HUNGARY TEL: +36 (1) 347 6070

E-MAIL: cyclolab@cyclolab.hu WEB: www.cyclolab.hu VAT No.: HU 10678970

ISO 9001

Fiber drums with locking clamp (Large, 100 L)

Material: fiber; Clamp material: Steel

Dimensions: 40 x 82.5 cm Total weight: 4.3 kg

Tests and methods	Requirements
Appearance (visual)	Free from visible mechanical contaminations and damages
Closure (visual)	Presence of sealing material

Storage: at room temperature.

Locking clamp (lead/seal): Mini Cable Lock Premium 1,5 mm, White with ADD with 30 cm wire, Metal or equivalent seal.

Cable tie: natural, white or black colored plastic minimum 100 mm long cable tie.

Desiccant: Active silica gel pad

Quality: meets the FDA requirements

Dimensions: 110 x 110 mm, weight: 100 g, 25 g

Tests and methods	Requirements
LOD (CoA of the manufacturer)	NMT 4.0 %
Water absorption capacity (RH 50%) (CoA of the manufacturer)	≥ 25.0 %

Storage: at room temperature, in tightly closed container

Packaging instructions: The dried SBECD bulk powder must be filled into a gamma sterilized PE bag that is twist-tied and quick shut with cable tie as the last step of the manufacturing. This PE bag is placed into the Triplex foil bag (secondary packaging material) with silica gel desiccant that is then heat sealed and labelled.

The packaged product is then transferred to HDPE containers or fiber drums as tertiary packaging material, and then secured with tamper evident seals.

Adequate labels should be placed on the inner PE bag, Triplex foil bag and HDPE drums for identification.



MAIL: BUDAPEST, P.O.BOX 435, H-1525 HUNGARY LOCATION: ILLATOS ÚT 7. BUDAPEST, H-1097 HUNGARY TEL: +36 (1) 347 6070

E-MAIL: cyclolab@cyclolab.hu

WEB: www.cyclolab.hu VAT No.: HU 10678970

ISO 9001



Blank packaging label (small):

CYCLO

CYCLOLAB MAR: BUDAPEST, P.O.BOX 435, H-1525 HUNGARY LOCATION: BLATOS UT 7, BUDAPEST, H-1097 HUNGARY TEL: -56 (f) 347 6070 E-MAIL: cyclolab (reyclolab bu WES: www.cyclolab bu VAT No.: HU 10678970

Product name:

Sulfobutyl-ether-\(\beta\)-cyclodextrin

sodium salt (SBECD)

Quality:

USP NF and EP (for manufacturing

of injectable dosage forms)

Batch No .:

pack

Net weight:

...... kg

Preserve in well-closed containers, store at room temperature.

Protect from moisture.

Expiry date:

Blank packaging label (big):

CYCLO PIONEERS IN CYCLODEXTRIN-BASED

CYCLOLAB

CYCLODEXTRIN RESEARCH & DEVELOPMENT LABORATORY LTD.

MAIL: BUDAPEST, P.O.BOX 435, H-1525 HUNGARY
LOCATION: ILLATOS UT 7. BUDAPEST, H-1097 HUNGARY
TEL: +36 (1) 347 6070

E-MAIL: cyclolab@cyclolab.hu WEB: www.cyclolab.hu VAT No.: HU 10678970

Product name:

Sulfobutyl-ether-β-cyclodextrin sodium salt

drum

Batch No .:

PO-No.:

Gross weight:

kg

Net weight:

..... kg

Storage:

Preserve in well-closed containers, store at room temperature.

Protect from moisture.

Expiry date:



MAIL: BUDAPEST, P.O.BOX 435, H-1525 HUNGARY LOCATION: ILLATOS ÚT 7. BUDAPEST, H-1097 HUNGARY TEL: +36 (1) 347 6070

E-MAIL: cyclolab@cyclolab.hu WEB: www.cyclolab.hu VAT No.: HU 10678970

ISO 9001



Photos:





£ 25 September, 2023 QC/date QA/date

Change

v02: "complies with EP and FDA regulation" was deleted, new tertiary packaging option: fiber drum, secondary packaging was modified with pressure relief

v03: one PE bag option, modified PE bag specification (weight criteria was deleted)

v04: the pictures were changed

v05: supplemented with cable tie and locking clamp

v06: updated the specification, corrected expiry date, addition of dimensions and total weights