

PRODUCT SPECIFICATION	
Product: 3,6-Anhydro-heptakis(6-deoxy-6-(2-carboxyethyl)thio)-gamma-cyclodextrin partly sodium salt (Mono-AH-SGM) Molecular formula: $C_{69}H_{106-n}Na_nO_{46}S_7$ Formula weight: $1896 + n \cdot 22$ Quality: analytical standard	Version: 01 Code: Rel_SGM_CY_3240_v01

Prepared by (QC)/date:	Revised by (QA)/date:	Approved by (QA)/date:
<i>Old Fozsai</i> 19 October 2022	<i>Mich. Póron</i> 19 October, 2022	<i>Old Fozsai</i> 19 October 2022

Tests and methods	Requirements
Appearance# / visual	white or off-white powder
Identification 1 / LC-MS (EP 2.2.29., 2.2.43)	conforms to theoretical mass
Identification 2 / NMR (EP 2.2.33)	conforms to theoretical structure
Identification 3 / IR (EP 2.2.24)	conforms to standard
Purity# ¹ / HPLC (EP 2.2.29)	> 90.0 area%
Other cyclodextrin-related impurities# ¹ / HPLC (EP 2.2.29)	< 10.0 area %
Sodium content ² / CE (EP 2.2.47)	< 3.0 w/w %
Residual solvents / TGM-MS (EP 2.2.34, 2.2.43)	< 10.0 w/w %
Purity calculated ³	to be reported

Tests to extend expiry

1 HPLC analysis to be carried out on an Agilent Poroshell column with H_3PO_4 / NH_3 buffer system at pH=2.6 and UV-Vis detection

2 The compound is supplied in slightly acidic form for stability reasons; the actual molar sodium content (n) is indicated on the Certificate of Analysis of each batch.

3 Calculated purity (w/w%) = ((100 (w/w%) - Sodium formate content (w/w%) - Residual solvents (w/w%)) * Purity (area%)) / 100

Notes: the compound is a single isomer.



CYCLOLAB
CYCLODEXTRIN RESEARCH & DEVELOPMENT
LABORATORY LTD.
MAIL: BUDAPEST, P.O.BOX 435, H-1525 HUNGARY
LOCATION: ILLATOS ÚT 7. BUDAPEST, H-1097 HUNGARY
TEL: +36 (1) 347 6070
E-MAIL: cyclolab@cyclolab.hu
WEB: www.cyclolab.hu
VAT No.: HU 10678970



Comments: only for reagent and laboratory purposes, not for human use.

Storage: in deep freezer, under inert atmosphere, in tightly closed container.

Retest date: 12 months from date of qualification.

<i>Version No:</i>	<i>Short description of change</i>
01	First version